

HOW TO APPLY THE INVISIBLE WALL HANGER TO YOUR BONE CHINA PLATE

- 1 First, clean the back of your new China Plate. Allow to dry completely before proceeding.
- 2 Using your fingers, add a couple drops of water to the glue on the back of the adhesive wall hanger. Leave for 5 to 10 minutes until the glue is tacky.
- 3 Apply the adhesive side of your wall hanger to the back of the plate. Ensure that the ring is in the correct position in the centre of the plate. Press firmly and rub well to get rid of any trapped air. Leave overnight to dry.
- 4 Before hanging, test the strength of the adhesive by pulling firmly on the ring.

To remove the wall hanger, soak the plate in warm water.

The invisible wall hanger is designed so that it can be used and removed safely without spoiling the plate or its markings.